

3-in-1 Weighing Scale



*Unique Digital Weighing Solution
for Infant, Child & Adult*

Infant weigh bed
up to 10 Kgs

Folding child seat
up to 40 Kgs

Adult standing
up to 120 Kgs

Digital Monitoring Device



Front View



Back View

* Battery backup for 6 hours



**Infant Digital Weighing Scale
Model IWS-101**



**3-in-1 Digital Weighing Scale
Model MWS-301**

* Technical specifications may be altered by Neokraft without prior notice, accessories shown in the catalogue are not part of the standard equipment

Manufactured by:



NEOKRAFT MEDICAL PRIVATE LIMITED

18, 1st Cross, 11th Main, Tata Nagar, Sahakara Nagar [P],
BANGALORE - 560 092, INDIA.

Phone : (080) 6500 0249

Mobile : +91 9686 78 5555

Email : info@neokraftmedical.com

Website : www.neokraftmedical.com

Channel Partner